



Magnetic Cell Selection and Separation of Human CD19+ Cells

OVERVIEW

The COL-iso™ Human PBMC CD19+ Cells Isolation Kit (Cat. # K10102) is designed to isolate CD19+ human PBMC cells using positive selection. The resulting cell preparation is highly enriched for CD19+ cells. Purity of recovered CD19+ cells can be up to 97%-99% and will vary depending on the preparation.



MATERIALS REQUIRED

1. Magnetic Separator e.g. MiniMacs Separator (Miltenyi Catalog # 130-042-102)
2. Column e.g. MS Column (Miltenyi, Catalog # 130-042-201)
3. Sterile serological and Pasteur pipettes or transfer pipettes
4. 30µM Filter (Partec, Catalog #04-0042-2316)
5. Bench top centrifuge
6. 2-8° C refrigerator
7. Deionized or distilled water

Cell Selection Principle

1. Positive selection of CD19+ cells is achieved by incubation with biotinylated anti-Human CD19 monoclonal antibody. CD19 monoclonal antibody bound cells are then
2. magnetically tagged with COL-iso™-Streptavidin.
3. Magnetically tagged CD19+cells are then retained in the magnetic column. (These are the desired cells); unwanted/untagged cells run through.
4. Upon removal of column from magnetic field, CD19+cells can then be eluted.

Additional Products and Services:

-  [Mouse Monoclonal Antibody](#)
-  [Rat Monoclonal Antibody](#)
-  [Rabbit Monoclonal Antibody](#)
-  [Human Monoclonal Antibody](#)
-  [Polyclonal Antibody](#)
-  [Antibody Sequencing](#)
-  [Hybridoma Sequencing](#)
-  [CAR T-cells](#)
-  [Lentivirus production](#)
-  [Cancer Stem Cells](#)
-  [Specialty Cell Culture Media](#)
-  [T-cell Expansion beads](#)

Ask about our full line of CRO services to provide supplemental assistance or the entire support necessary to complete your project on time and with the data you need to move forward.

Magnetic Cell Selection and Separation of Human CD19+ Cells



Cell Selection Capacity

Separator	Max No. of CD19 ⁺ cells/column	Max No. of cells/column
Column/EA	*1x10 ⁷	*1x10 ⁸

*: The Max No. of cells will vary by $\pm 40\%$ depending on the preparation

1. Biotinylated anti-Human CD19 Antibody (Part C10102) - 2mL
2. COL-isoTM -Streptavidin. (Part B10002) -2mL proprietary formulation.
3. DRNase (proprietary formulation of DNase I and RNase) – 1mL (Part DR10100)

O

k $\frac{1}{\sqrt{2}}$ kV
 k $\frac{1}{\sqrt{2}}$ kV
 o $\frac{1}{\sqrt{2}}$ kV
 o $\frac{1}{\sqrt{2}}$ kV

Cell Selection Procedure

- I. **Cell Preparation:** Cells and reagents should be kept cold using an ice bath or a refrigerator unless otherwise specified. Incubations must be carried out at 2-8°C in a refrigerator and not in an ice bath to avoid excessively low temperatures that can slow the kinetics of the optimized reactions.
- A. **Preparing a cell suspension from frozen PBMC/mPBMC/CBMNC**
 1. To a 50 mL conical tube add 10µL formulated DRNase per 10⁷ cells.
 2. Transfer desired amount of cell suspension to the 50ml conical tube.
 3. Drop wise add 15mL pre-warmed (37°C) DMEM containing 10% FBS to the cells with constant swirling.
 4. Centrifuge cell suspension at 300 x g at 4°C for 15 minutes.
 5. Carefully remove all but approximately 100µL of the supernatant using a pipette.
 6. Gently resuspend 10⁷cells with 80uL COLD Buffer.
 7. Pre - wet a 30-50µm nylon cell strainer then pass the suspended cells through the strainer.



Magnetic Cell Selection and Separation of Human CD19+ Cells

B. Preparing a cell suspension from fresh PBMC/mPBMC/CB MNC

1. Centrifuge cell suspension at 300 x g at 4°C for 15 minutes.
2. Gently resuspend 10^7 cells with 80 μ L COLD Buffer.
3. Pre-wet a 30-50 μ m nylon cell strainer then pass the suspended cells through the strainer.

Cells must be resuspended in cold reaction buffer prior to the antibody selection procedure. Buffer has to be kept on ice at all times.

NOTE: For downstream applications that are sensitive to DRNase (eg. hematopoietic colony assays), wash cells once in the appropriate assay buffer (without DRNase) before continuing.

II. Magnetic labelling of CD19+ cells

1. Transfer desired amount PBMC cells to an Eppendorf tube.
2. Add 20 μ L of biotinylated anti-human CD19 antibody (Part C10102) per 10^7 cells.
3. Gently mix the cell-antibody suspension, avoiding formation of bubbles, and incubate at 2-8° C on a rotator for 15 minutes.
4. After incubation, wash cells by adding 1-2 mL of buffer per 10^7 cells and centrifuge at 4°C at 300 x g for 10 minutes.
5. Carefully remove supernatant and resuspend 10^7 cells in 80 μ L of buffer.
6. Add 20 μ L COL-iso™-Streptavidin (Part B10002) per 10^7 cells.
7. Mix gently and incubate at 2- 8° C on a rotator in a refrigerator for 15 minutes.
8. After incubation, wash cells by adding 1-2 mL of buffer per 10^7 cells and centrifuge at 4°C at 300 x g for 10 minutes.
9. Completely remove supernatant and gently resuspend cell pellet up to 10^8 cells in 500 μ L of buffer.

III. Magnetic Separation

- IV. Place column in magnetic field. Prime column by rinsing 1x with 500 μ L of filtered Buffer.
- V. Load up to 10^8 cell suspension onto each equilibrated column. (i.e. 2×10^8 cells would require the use of 2 MS columns) Carefully save effluent as Flow Through.
- VI. Wash column 3x with 500 μ L of cold buffer. Only apply new buffer when column reservoir is empty. Collect effluent into Flow Through from step 2.





Magnetic Cell Selection and Separation of Human CD19+ Cells



4. At the end of the washing step, remove column from magnetic field and place column on a collection tube.
5. Add 1mL of buffer onto column and immediately flush out the CD19+cells with plunger. Label tube as Elution.
6. Centrifuge Flow Through and Elution at 4°C at 300 x g for 5 minutes.
7. Cells are now ready for further experimentation or FACS analysis.